

TONO-PEN AVIA[®] VET[™] TONOMETER

USER MANUAL



Item no. 173600

TELEPHONE +49 7461 96 580 0

www.eickemeyer.com

veterinary technology for life
Eickemeyer[®]

The Tono-Pen AVIA Tonometer Help Line

Should you need immediate help with a technical question or guidance through the appropriate procedure, just call the Help Line at 888-849-8955

Note: When contacting our Technical Service Group at Reichert, please have the appropriate product number, product serial number, date of purchase, and nature of inquiry available.

Product Number: _____

Serial Number: _____

Date of Purchase: _____

©2017 AMETEK, Inc.

Reichert, Reichert Technologies, Tono-Pen, Tono-Pen AVIA, Tono-Pen AVIA POWERCEL and Ocu-Film are registered trademarks of Reichert, Inc.

AMETEK is a registered trademark of AMETEK, Inc.

All other trademarks are property of their respective owners.

The information contained in this document was accurate at time of publication. Specifications subject to change without notice. Reichert, Inc. reserves the right to make changes in the product described in this manual without notice and without incorporating those changes in any products already sold.

ISO 9001/13485 Certified – Reichert products are designed and manufactured under quality processes meeting ISO 9001/13485 requirements.

No part of this publication may be reproduced, stored in a retrieval system, or transmitted in any form or by any means, electronic, mechanical, recording, or otherwise, without the prior written permission of Reichert, Inc.

Caution: Federal law restricts this device to sale by or on the order of a licensed physician. Rx only.

Table of Contents

Warnings and Cautions	4
Symbols.....	8
Introduction.....	9
Indications for Use	9
Contraindications	9
Device Description	10
Features	10
Parts Identification.....	11
Available Accessories	11
Instructions for Use	12
Tono-Pen AVIA Tonometer Modes	12
Tono-Pen AVIA Tonometer Verification.....	13
Tono-Pen AVIA Tonometer Tones.....	15
Battery Installation and Replacement	16
Liquid Crystal Display - User Interface.....	17
Tono-Pen AVIA Tonometer Preparation.....	19
Verification Mode.....	20
Performing IOP Measurements.....	21
Applanation	22
Interpreting the IOP Measurement.....	23
Cleaning & Maintenance	24
Cleaning Instructions.....	24
Battery	25
Storage.....	25
Disposal	25
Troubleshooting.....	26
Specifications	27
Guidance & Manufacturer’s Declaration	28
Device Regulatory Classification.....	30
Warranty	31

Warnings and Cautions

Reichert Technologies (Reichert) is not responsible for the safety and reliability of this instrument when:

- Assembly, disassembly, repair or modification is made by unauthorized dealers or persons.
- Instrument is not used in accordance with its User's Guide.



WARNING: AN INSTRUCTION THAT DRAWS ATTENTION TO RISK OF INJURY OR DEATH.

WARNING: UNITED STATES FEDERAL LAW AND EUROPEAN REGULATIONS REQUIRE THAT THIS DEVICE BE PURCHASED ONLY BY A PHYSICIAN OR A PERSON ACTING ON BEHALF OF A PHYSICIAN.

WARNING: DO NOT REPAIR OR SERVICE THIS INSTRUMENT WITHOUT AUTHORIZATION FROM THE MANUFACTURER. ANY REPAIR OR SERVICE TO THIS INSTRUMENT MUST BE PERFORMED BY EXPERIENCED PERSONNEL OR DEALERS WHO ARE TRAINED BY REICHERT OR SERIOUS INJURY TO THE OPERATOR OR PATIENT MAY OCCUR.

WARNING: THIS INSTRUMENT SHOULD BE USED IN STRICT ACCORDANCE WITH THE INSTRUCTIONS OUTLINED IN THIS USER'S GUIDE. THE SAFETY OF THE OPERATOR AND THE PERFORMANCE OF THE INSTRUMENT CANNOT BE GUARANTEED IF USED IN A MANNER NOT SPECIFIED BY REICHERT TECHNOLOGIES.

WARNING: MODIFICATIONS TO THIS INSTRUMENT ARE NOT ALLOWED. ANY MODIFICATION TO THIS UNIT MUST BE AUTHORIZED BY REICHERT OR SERIOUS INJURY TO THE OPERATOR OR PATIENT MAY OCCUR.

WARNING: IF THIS INSTRUMENT IS MODIFIED, APPROPRIATE INSPECTION AND TESTING MUST BE CONDUCTED TO ENSURE CONTINUED SAFE USE OF THIS INSTRUMENT.

WARNING: DO NOT USE THE TONO-PEN AVIA TONOMETER ON A PATIENT WITHOUT AN OCU-FILM + TIP COVER OR INACCURATE READINGS MAY BE OBTAINED.

WARNING: DO NOT USE AN OCU-FILM + TIP COVER ON MORE THAN ONE PATIENT TO HELP PREVENT CROSS CONTAMINATION.

WARNING: IT IS IMPERATIVE THAT A FRESH OCU-FILM + TIP COVER BE USED FOR EACH PATIENT FOR THE PROTECTION OF THE EYE AND THE TRANSDUCER ASSEMBLY. THE TRANSFER OF INFECTION IS THUS PREVENTED. REPLACE OCU-FILM + TIP COVERS AFTER EXAMINATION OF EACH PATIENT.

Warnings and Cautions (continued)

WARNING: DO NOT USE THE TONO-PEN AVIA TONOMETER IF THE TRANSDUCER ASSEMBLY IS CRACKED, CHIPPED OR SHOWS ANY IRREGULARITY OF THE SURFACE, TO PREVENT PATIENT INJURY, AND/OR INACCURATE READINGS.

WARNING: DO NOT USE EXCESSIVE PRESSURE DURING APPLANATION OR EYE INJURY MAY OCCUR.

WARNING: THE CORNEAL SURFACE NEEDS ONLY TO BE MOMENTARILY TAPPED, INDENTATION OR ADDITIONAL PRESSURE AFTER THE AUDIO “CHIRP” TONE IS HEARD IS NOT REQUIRED AND MAY LEAD TO DAMAGE TO THE EYE. IF MEASUREMENT PROCESS CAUSES PATIENT DISCOMFORT, DEVICE DOES NOT SHOW A MEASUREMENT, OR IF DEVICE REQUIRES MULTIPLE ATTEMPTS TO OBTAIN A MEASUREMENT, STOP THE EXAMINATION AND REFER TO THE TROUBLESHOOTING SECTION OF THIS MANUAL.

WARNING: OCU-FILM + TIP COVERS CONTAIN NATURAL LATEX WHICH MAY CAUSE ALLERGIC REACTIONS. QUESTION PATIENTS ABOUT ALLERGIES TO LATEX BEFORE MEASURING THEM WITH THE TONO-PEN AVIA TONOMETER.

WARNING: DO NOT CARRY THE TONO-PEN AVIA POWERCEL BATTERY IN A POCKET, OR CLOSE TO YOUR PERSON, AS A BURN INJURY MAY RESULT.

WARNING: THE BATTERY SHOULD ONLY BE REPLACED WITH THE BATTERY SPECIFIED IN THIS MANUAL. USE OF ANOTHER BATTERY MAY CAUSE FIRE OR AN EXPLOSION.

WARNING: DO NOT PLACE A SHORTING DEVICE BETWEEN THE BATTERY TERMINALS, OR ALLOW THE BATTERY TO BECOME WET. MISUSE OR IMPROPER DISPOSAL OF THIS BATTERY MAY CAUSE IT TO BECOME VERY HOT, IGNITE OR EXPLODE. DAMAGE TO THIS UNIT AND/OR SERIOUS PERSONAL INJURY MAY RESULT.

WARNING: DO NOT RECHARGE THE BATTERY. THE BATTERY IS NOT DESIGNED TO BE CHARGED BY ANY ELECTRICAL SOURCE. CHARGING COULD GENERATE GAS AND INTERNAL SHORT-CIRCUITING, LEADING TO DISTORTION, LEAKAGE, OVERHEATING, EXPLOSION OR FIRE.

WARNING: DO NOT EXPOSE THE BATTERY TO TEMPERATURES ABOVE 140°F, DISASSEMBLE THE BATTERIES, OR DAMAGE TO THIS UNIT AND/OR SERIOUS PERSONAL INJURY MAY RESULT.

-continued-

Warnings and Cautions (continued)

WARNING: NEVER ALLOW LIQUID LEAKING FROM THE BATTERY TO GET IN YOUR EYES OR MOUTH AS THIS LIQUID COULD CAUSE SERIOUS PERSONAL INJURY. IF IT COMES IN CONTACT WITH YOUR EYES OR MOUTH, FLUSH THEM IMMEDIATELY WITH PLENTY OF WATER AND CONSULT A PHYSICIAN.

WARNING: ALWAYS KEEP BATTERIES OUT OF THE REACH OF INFANTS AND YOUNG CHILDREN TO PREVENT THEM FROM BEING SWALLOWED. IF SWALLOWED, CONSULT A PHYSICIAN IMMEDIATELY.

WARNING: THE USE OF ACCESSORIES OR CABLES OTHER THAN THOSE SPECIFIED, WITH THE EXCEPTION OF THOSE SOLD BY THE MANUFACTURER AS REPLACEMENT PARTS FOR INTERNAL COMPONENTS, MAY RESULT IN INCREASED EMISSIONS OR DECREASED IMMUNITY OF THE EQUIPMENT OR SYSTEM.



CAUTION: AN INSTRUCTION THAT DRAWS ATTENTION TO THE RISK OF DAMAGE TO THE PRODUCT.

CAUTION: DO NOT USE THE TONO-PEN AVIA TONOMETER WITHOUT AN OCU-FILM + TIP COVER TO PREVENT DAMAGE TO THE TRANSDUCER ASSEMBLY.

CAUTION: DO NOT TOUCH THE TRANSDUCER ASSEMBLY WITHOUT AN OCU-FILM + TIP COVER APPLIED OR DAMAGE TO THE TRANSDUCER ASSEMBLY MAY RESULT.

CAUTION: DO NOT BUMP, JAR OR DROP THE DEVICE, OR DAMAGE TO THE ELECTRONICS MAY OCCUR.

CAUTION: DO NOT IMMERSE THE TONO-PEN AVIA TONOMETER IN FLUIDS OR DAMAGE TO THE ELECTRONICS MAY OCCUR.

CAUTION: DO NOT USE SOLVENTS OR STRONG CLEANING SOLUTIONS ON ANY PART OF THIS INSTRUMENT AS DAMAGE TO THE UNIT MAY OCCUR. SEE MAINTENANCE SECTION FOR DETAILED CLEANING INSTRUCTION.

CAUTION: USE OF AMMONIA BASED CLEANERS ON THE LIQUID CRYSTAL DISPLAY (LCD) MAY CAUSE DAMAGE TO THE DISPLAY. SEE MAINTENANCE SECTION FOR DETAILED CLEANING INSTRUCTION.

CAUTION: DO NOT ATTEMPT TO STERILIZE THE TONO-PEN AVIA TONOMETER OR DAMAGE TO THE ELECTRONICS MAY OCCUR.

Warnings and Cautions (continued)

CAUTION: DO NOT AUTOCLAVE OR DISINFECT USING HIGH TEMPERATURES EXCEEDING THE RECOMMENDED TEMPERATURES INDICATED IN THE SPECIFICATIONS SECTION OF THIS MANUAL OR DAMAGE TO THE UNIT MAY OCCUR.

CAUTION: DO NOT ATTEMPT TO MODIFY THE TONO-PEN AVIA TONOMETER OR THE TONO-PEN AVIA POWERCEL BATTERY OR DAMAGE TO THE DEVICE MAY OCCUR.

CAUTION: DO NOT STORE THE TONO-PEN AVIA TONOMETER WITHOUT AN OCU-FILM + TIP COVER OR DEBRIS MAY ENTER THE TRANSDUCER ASSEMBLY AND CAUSE MALFUNCTIONS.

CAUTION: MEDICAL ELECTRICAL EQUIPMENT NEEDS SPECIAL PRECAUTIONS REGARDING EMC AND NEEDS TO BE INSTALLED AND PUT INTO SERVICE ACCORDING TO THE EMC INFORMATION PROVIDED IN THIS GUIDE. PORTABLE AND MOBILE RF COMMUNICATIONS EQUIPMENT CAN AFFECT MEDICAL ELECTRICAL EQUIPMENT.

CAUTION: ELECTROMAGNETIC INTERFERENCE FROM OTHER DEVICES MAY AFFECT THIS INSTRUMENT. IF INTERFERENCE IS PRESENT, TURN OFF OTHER ELECTRONIC DEVICES, OR REMOVE THEM FROM THE IMMEDIATE AREA WHILE OPERATING THIS INSTRUMENT.

CAUTION: PORTABLE AND MOBILE RF COMMUNICATIONS EQUIPMENT CAN EFFECT MEDICAL ELECTRICAL EQUIPMENT.

CAUTION: THIS INSTRUMENT IS NOT TO BE USED NEAR HIGH-FREQUENCY EMITTING SURGICAL EQUIPMENT.

CAUTION: OCU-FILM + TIP COVERS SHOULD BE STORED BETWEEN 35° AND 80° FAHRENHEIT (2°-27° CELSIUS).

CAUTION: EACH BOX OF HIGH QUALITY OCU-FILM + TIP COVERS HAS A "USE BEFORE" DATE STAMPED ON THE BOX. WE SUGGEST USE OF OUR OCU-FILM + BEFORE THIS DATE TO GUARD AGAINST THE POSSIBILITY OF SHELF WEAR. THE LATEX MATERIAL USED IN THE OCU-FILM + TIP COVERS CAN DEGRADE. A DEGRADED FILM MAY RESULT IN LEAKAGE OF WETTING SOLUTION AND INSTRUMENT DAMAGE. EXAMINE EACH OCU-FILM + TIP COVER FOR YELLOWING, CRACKS, OR A STICKY TEXTURE PRIOR TO USE.

Symbols



Caution

REF

Catalog Number

SN

Serial Number



Date of Manufacture



Manufacturer



Waste of Electrical and Electronic Equipment



Compliance to Medical Device Directive 93/42/EEC



Consult Instructions for Use



Authorized Representative in European Community



Fragile Contents in Shipping Container - handle with care



Do not get Shipping Container wet



Type BF Applied Part

SYMBOLS FOR OCU-FILM + TIP COVERS ONLY



Do not reuse. Single Use



Use By

LOT

Lot Number



Contains natural rubber latex

Introduction

Congratulations on your purchase of the Tono-Pen AVIA® tonometer.

The Tono-Pen AVIA tonometer is a prescription only device intended for measuring intraocular pressure (IOP) during routine eye examinations or when increased intraocular pressure is suspected by properly trained eyecare professionals such as ophthalmologists, optometrists, opticians, and eye care technicians.

This User's Guide is designed as a training and reference manual for operation, maintenance, and troubleshooting. We recommend that you read it carefully prior to use and follow the instructions in the guide to ensure optimum performance of your new instrument. If used properly, the Tono-Pen AVIA tonometer will provide you with fast, accurate and reliable measurements for many years. Properly trained eyecare professionals such as ophthalmologists, optometrists, opticians and eye care technicians should operate this instrument.

Please retain this manual for future reference and to share with other users. For additional copies of this manual or questions related to the Tono-Pen AVIA tonometer, contact your local authorized Reichert® dealer or contact our Customer Service department directly at:

Tel: 716-686-4500

Toll Free: 888-849-8955

Fax: 716-686-4555

E-mail: reichert.information@ametek.com

Indications for Use

The indications for use include measuring intraocular pressure (IOP) for suspected glaucoma, or when increased intraocular pressure is suspected.

Contraindications

None.

Device Description

The Tono-Pen AVIA tonometer is an ergonomic, hand-held tonometer that measures intraocular pressure. The body of the instrument is designed to fit comfortably in the user's hand, facilitating fast and accurate measurements. The tip of the Transducer Assembly converts applied force into an electrical signal. The electronics housed in the ergonomic Tono-Pen AVIA tonometer body process and analyze the waveforms produced by each applanation of the corneal surface of the eye. These are used to produce an averaged IOP measurement. The measurement is displayed on the Liquid Crystal Displays (LCDs).

A replaceable battery compartment houses the Tono-Pen AVIA POWERCEL® Battery, consisting of two Lithium Manganese Dioxide batteries.

Features

The Tono-Pen AVIA tonometer has the following features:

- Easy to use - IOP can be measured reliably by medical professionals.
- Portable - The Tono-Pen AVIA tonometer weighs just 71 g (2.4 oz) and is battery operated.
- Accurate - The measurements from the Tono-Pen AVIA tonometer correlate strongly with Goldmann applanation tonometry and direct measurements of IOP.
- Versatile - The Tono-Pen AVIA tonometer may be used easily with the patient in any position, making the instrument suitable for the office, in glaucoma clinics, at the hospital bedside, and in remote locations.

Introduction (continued)

Parts Identification

- A** Tip - tonometer tip
- B** Transducer Assembly - housing for tonometer tip
- C** LCD - displays the IOP in mmHg, number of applanations collected, statistical confidence indicator, and battery life status
- D** LED - green light on in Applanation Mode.
- E** Tono-Pen AVIA POWERCEL Battery - battery
- F** Activation Button - Applanation Mode select button.



Figure 1. Key Features

Accessories

- | | |
|---------|--------------------------------------|
| 230577 | Tono-Pen AVIA POWERCEL Battery |
| 230651 | Ocu-Film® + Tip Covers (150 per box) |
| 68C1334 | Carrying Case |
| 68E3892 | User's Guide |
| 980051 | Quick Reference Guide |

Instruction for Use

Tono-Pen AVIA Tonometer Modes

Power Up Mode

Pressing the activation button initiates the battery life, LCD, Transducer Assembly, and electronic self-test. A successful power-up puts the device automatically into applanation mode.

Applanation Mode

A 15 second interval within which the user should complete the 10 applanations needed to acquire an IOP value. During applanation mode, the LCD will display the number of applanations collected. After 10 applanations are collected, the LCD will display the IOP along with a statistical confidence indicator. After 25 seconds, the device will automatically go into a sleep mode. If the activation button is pressed again the data will be lost and the unit will be ready to start another measurement sequence.

Verification Mode

Pressing and holding the activation button for 5 seconds initiates an electronic self test mode. See [Tono-Pen AVIA Tonometer Verification](#) section for details.

Sleep Mode

A power-saving mode automatically activated after 25 seconds of non-use. Pressing the activation button will initiate a power-up mode.

Instruction for Use (continued)

Tono-Pen AVIA Tonometer Verification

If suspect readings are observed, a verification test should be performed. This will ensure that the Transducer Assembly and electronics are performing correctly.

Verification

1. Hold the Tono-Pen AVIA tonometer with the Transducer Assembly end pointing down towards the floor.
2. Press and hold the activation button for 5 seconds - a beep will sound at one second intervals.
3. At the end of the 5 second button hold, the display will show **[dn]**.
4. Keep the pen vertical, with the Transducer Assembly pointing down towards the floor, for a total of 15 seconds.
5. At the end of this period, a beep will sound and the display will show **[UP]**.
6. Immediately point the Transducer Assembly straight up and wait for the next beep (within 3 seconds).



Figure 2. Tono-Pen AVIA Verification Procedure
-continued-

Instruction for Use (continued)

Verification (continued)

7. A properly functioning Tono-Pen AVIA tonometer will display **[Pass]**. Pressing the activation button will now put the device into Applanation Mode.



Figure 3. Tono-Pen AVIA Verification Pass

8. Verification test will need to be repeated if the display shows **[FAIL]**. A single press of the Activation button will re-start the verification test.

Note: The Tono-Pen AVIA tonometer will remain in the Verification mode until it passes verification test.

9. The most common cause for failure is debris in the tip of the Transducer Assembly. Follow the instructions for cleaning the tip in the “Maintenance & Storage” section before re-performing the verification test.
10. If this is unsuccessful, please contact the Reichert Technical Service Department.

Instruction for Use (continued)

Tono-Pen AVIA Tonometer Tones

The Tono-Pen AVIA tonometer generates two different tones.

A “BEEP” will sound when:

- Power Up Mode fails
- Verification Mode fails
- Applanation Mode starts
- IOP testing is completed
- IOP values are out of range (<5 or >55)
- At various points during the verification test.

A “CHIRP” will sound when:

- A valid IOP measurement has been taken during applanation.

Instruction for Use (continued)

Battery Installation and Replacement

The Tono-Pen AVIA tonometer is supplied with a Tono-Pen AVIA POWERCEL Battery that needs to be installed prior to use. The Tono-Pen AVIA POWERCEL Battery is the only replacement battery that can be used with this device.


1. Insert the Tono-Pen AVIA POWERCEL Battery into the base of Tono-Pen AVIA tonometer until fully seated. The Battery is “keyed” for correct installation.
2. Press the Activation Button.
3. Check that the battery symbol in the LCD has all segments displayed.
4. Check that the Self-Test is initiated.
5. Dispose of used Batteries in accordance with local regulations.





Figure 4. Tono-Pen AVIA POWERCEL Battery Installation and Replacement




Instruction for Use (continued)

Liquid Crystal Display - User Interface

POWER UP MODE - USER INTERFACE SEQUENCE OF EVENTS		
	User Interface LCD/Battery Test	LED Off Passed Battery Life Test/Passed LCD Test




(SEE APPLANATION MODE)




POWER UP MODE - USER INTERFACE ERROR CODES		
	User Interface Error Code	LED Off Battery Level Critical (when low, symbol will flash)
	User Interface Error Code	LED Off Failed Self Test


APPLANATION MODE - USER INTERFACE SEQUENCE OF EVENTS		
	User Interface Ready to Measure	LED On, 1 Beep Ready to Begin IOP Testing (15 seconds for testing before Time Out)
	User Interface Testing Display	LED On During Test, screen shows number of applanations achieved (in this case 4 of 10) Each applanation taken equals 1 Chirp
	User Interface Test Complete	LED Off, 1 Beep Test Complete Patient has an IOP of 16 with a statistical confidence indicator of 95

Instruction for Use (continued)

Liquid Crystal Display - User Interface (continued)

APPLANATION MODE - USER INTERFACE ERROR CODES		
	User Interface Error Code	LED Off Error - IOP value Under Range (<5mm Hg)
	User Interface Error Code	LED Off Error- IOP Value Over Range (>55mm Hg)
	User Interface Error Code	LED Off Timeout - not enough IOP values collected

VERIFICATION MODE - SEQUENCE OF EVENTS		
	User Interface Begin Verification Mode	Hold Tono-Pen AVIA tonometer with transducer down towards floor (LED Off)
	User Interface Invert Pen	Invert Tono-Pen AVIA tonometer (LED Off)
	User Interface Successful Verification	Successful Tono-Pen AVIA tonometer verification (LED Off)

VERIFICATION MODE - USER INTERFACE ERROR CODE		
	User Interface Error Code	Unsuccessful Tono-Pen AVIA tonometer verification (LED Off)

Instruction for Use (continued)

Tono-Pen AVIA Tonometer Preparation

1. Allow the instrument to thermally stabilize to room temperature for approximately 30 minutes prior to use.
2. Remove the storage Ocu-Film + tip cover from the Transducer Assembly.
3. Visually inspect the Transducer Assembly for cracks, chips or other irregularities. Do not use if these conditions are present.
4. Slide an Ocu-Film + Tip Cover onto the Transducer Assembly until the ridge is seated, taking care not to apply the tip cover too tightly or too loosely. Make certain that the rubber is flat across the tip, but not taut.



Figure 5. Tono-Pen AVIA Tonometer Preparation

Note: Use of the Ocu-Film + Tip Cover is required under original Tono-Pen instrument warranty and service contract terms.

Note: Corneal topical anesthetic is required with tonometry.

Instruction for Use (continued)

Verification Mode

A verification Test may be performed to ensure that the Transducer Assembly and electronics are performing as expected. A verification test may also be initiated if suspect readings are observed. (See Tono-Pen AVIA Tonometer Verification section for details.)

If the verification test resulted with the LCD displaying **[FAIL]**, repeat the verification test, or reference the troubleshooting section for potential causes for test failure.

Note: The Tono-Pen AVIA tonometer will remain in the Verification mode until it passes verification.

Instruction for Use (continued)

Performing IOP Measurements

WARNING: OCU-FILM + TIP COVERS CONTAIN NATURAL RUBBER LATEX WHICH MAY CAUSE ALLERGIC REACTIONS. QUESTION PATIENTS ABOUT ALLERGIES TO LATEX BEFORE MEASURING THEM WITH THE TONO-PEN AVIA TONOMETER.

WARNING: DO NOT USE AN OCU-FILM + TIP COVER ON MORE THAN ONE PATIENT TO HELP PREVENT CROSS CONTAMINATION.

Note: Refer to the Reichert website for instructional videos.

1. Instill a drop of topical anesthetic into the eye to be examined.
2. Position the patient, seated or supine, in front of a fixation target.

Note: The Tono-Pen AVIA tonometer will function in any orientation.

3. Instruct the patient to look straight ahead at the fixation target (i.e. ear, nose, distant object) to minimize eye movement, with eyes fully open.
4. Hold the Tono-Pen AVIA tonometer as you would a pencil and position yourself to allow viewing of the Transducer Assembly and the patient's cornea where contact will be made. For normal corneas, central corneal contact is recommended.
5. The corneal surface needs only to be momentarily contacted. **Indentation or excessive force is not required and may lead to inaccurate readings or patient injury.** Refer to Figure 6.



Figure 6. Corneal Applanation

6. Brace the heel of your hand on the patient's cheek for stability while holding The Tono-Pen AVIA tonometer perpendicular to and within approximately 1.3 cm (1/2") of the patient's cornea.

Note: In the unlikely event an Electrostatic Discharge (ESD) occurs, the Tono-Pen AVIA tonometer may display inaccurate data. Simply re-press the Activation button and start the measurement again.

Instruction for Use (continued)

Applanation

1. Press and release the activation button once to initiate the applanation mode.
2. A brief display of [88888] will flash on the LCD, a double row of dashes [==] will be displayed, the green LED will turn on, and a “beep” tone will sound, indicating the Tono-Pen AVIA tonometer is ready to measure intraocular pressure.
3. Tap very lightly and briefly on the corneal surface. The device will chirp and the DATA field on the LCD will increment for each valid IOP reading obtained.
4. When another beep tone is heard, indicating 10 applanations have been read, the green LED will turn OFF, and the averaged IOP measurement will appear on the LCD above “mmHg”. The statistical confidence indicator will appear on the LCD above DATA.

Note: If at least 6 applanations were read, the IOP will be displayed after a 4 second delay, along with the statistical confidence indicator.

Note: After pressing the Activation button, do not shake the Tono-Pen AVIA tonometer as it may register a reading. If this occurs simply re-press the Activation button and start the measurement again.

If 15 seconds elapse prior to applanation being initiated, the Tono-Pen AVIA tonometer will initiate the Sleep Mode, indicated initially by the display of a single row of dashes [---], followed by a blank display. The applanation mode may be initiated by pressing the activation button.

If any error codes appear on the LCD after the final beep, the applanation procedure must be repeated.

Instruction for Use (continued)

Interpreting the IOP Measurement

During the measurement, the LCD shows the cumulative number of applanations detected. Once 10 applanations are achieved, the LCD will display the IOP in millimeters of mercury (mm Hg), along with a statistical confidence indicator.

A statistical confidence indicator of 95 means that the standard deviation of the valid measurements is 5% or less of the number shown. The higher the statistical confidence indicator, the more reliable the measurement.

If the statistical confidence indicator is 80 or 80-, a repeat measurement is recommended.



Figure 7. LCD showing an IOP measurement of 16 and statistical confidence indicator of 95

Cleaning & Maintenance

Cleaning Instructions

The Tono-Pen AVIA tonometer may have difficulty taking measurements or display **[FAIL]** after a verification when its tip is dirty and requires cleaning. When the Transducer Assembly of the tonometer has dirt and contaminants in the tip, cleaning of the tip is necessary. When the tip contains contaminants, it cannot move freely and the Tono-Pen AVIA may have erratic readings and then show a **[FAIL]** verification.



Figure 8. Cleaning the Tip

1. Remove Ocu-Film + tip cover from the tonometer, if one is installed.
2. Using canned air, place the tip of the Transducer Assembly against the outlet of the canned air as shown in Figure 8.
3. Blow the canned air into the tip of the tonometer for approximately 3 seconds.

Note: It is necessary to blow canned air directly into the tip so that the contaminants are pushed out.

4. After cleaning the tip with compressed air, the Transducer Assembly will be cold. Allow the Transducer Assembly to warm to room temperature.
5. Perform the tonometer verification as indicated in the [Tono-Pen AVIA Tonometer Verification](#) section of this manual.

Note: If the tonometer does not pass the verification procedure, then repeat the above cleaning instructions. Do not clean more than 3 times in a row. If the tonometer still will not pass verification, contact Reichert.

Note: Never use the Tono-Pen AVIA tonometer without an Ocu-Film + tip cover installed.

Note: Ocu-Film + tip covers are the only manufacturer approved covers for use with the Tono-Pen AVIA tonometer. Use of any other type of branded tip cover may affect readings or may void your warranty.

Cleaning & Maintenance (continued)

Cleaning Instructions (continued)

Note: Always store the Tono-Pen AVIA tonometer with an Ocu-Film + tip cover installed to protect the tonometer tip from dirt and contaminants.

Suggested Cleaning Schedule

Number of Patients per Week	Number of Days Between Cleaning
10	30
100	15
300	7
600	1

Battery

Replace the Tono-Pen AVIA POWERCEL Battery when:

- The battery indicator shows low battery after depressing the activation button.
- There are no beeps, the LCD remains blank after pressing the activation button, or a noticeable slowing occurs when activating the device.

Note: Replace the battery with P/N 230577.

Storage

- Cover the tip with an Ocu-Film + tip cover for protection.
- If the Tono-Pen AVIA tonometer is not to be used until a subsequent day, the instrument and accessories should be placed in the storage case provided.
- If the instrument is to be stored for an extended period, remove the Tono-Pen AVIA POWERCEL Battery to avoid possible damage to the instrument due to battery leakage.

Disposal


The Tono-Pen AVIA, Ocu-Film, and Tono-Pen AVIA POWERCEL do not generate any environmentally hazardous residues. At the end of its product service life, follow your local laws and ordinances regarding the proper disposal of this equipment.

Troubleshooting

The table below provides a guide for troubleshooting some basic Tono-Pen AVIA tonometer operational problems. If a problem persists after using this guide contact Reichert technical services.

SYMPTOM	PROBABLE CAUSE	CORRECTION
Battery symbol flashes	Low Tono-Pen AVIA POWERCEL Battery capacity	Replace Tono-Pen AVIA POWERCEL Battery
Multiple inaccurate readings	Improper technique (Example: too much or too little force used, taps too long, device not held perpendicular to the patient's eye)	Review "APPLANATION" Section of this manual
	Old or improperly applied Ocu-Film + tip cover being used, or non-Reichert tip cover being used	Replace Ocu-Film tip cover
	Debris in tip	Clean tip
	Mechanical or electronic damage	Arrange for service through Reichert Technical Service Group
No beeps and/ or no dashes upon activation	Activation button not properly pressed	Press Activation button
	No Tono-Pen AVIA POWERCEL Battery capacity	Replace Tono-Pen AVIA POWERCEL Battery
	Mechanical or electronic damage	Arrange for service through Reichert Technical Service Group
Verification Failure	Debris in tip	Clean tip
	Compressed air cleaning has lowered the temperature of the instrument	Allow Tono-Pen AVIA tonometer unit to warm to room temperature

Specifications


<p>PHYSICAL DIMENSIONS</p> <p>Size: 16 x 2 x 4.4 cm (6 1/4" x 3/4" x 1 3/4")</p> <p>Weight: 71 g (2.4 oz)</p>	<p>ENVIRONMENTAL REQUIREMENTS</p> <p>Operational Environment Ambient Temperature range: 15° to 35°C (59° to 95°F)</p> <p>Relative Humidity range: 30 to 75%</p> <p>Atmospheric Pressure range: 80 kPa to 106 kPa (23.6 to 31.3 in.Hg)</p> <p>Transport and Storage Environment Ambient Temperature range: -10° to 60°C (14° to 140°F)</p> <p>Relative Humidity range: 10 to 80% (non-condensing)</p> <p>Atmospheric Pressure range: 50 kPa to 106 kPa (14.8 to 31.3 in.Hg)</p>								
	<p>RANGE OF IOP MEASUREMENTS Accuracy by manometric measurements (95% Confidence)</p> <table border="1" data-bbox="477 1094 1003 1257"> <thead> <tr> <th>Measurement (mmHg)</th> <th>Accuracy (mmHg)</th> </tr> </thead> <tbody> <tr> <td>5-25</td> <td>± 1.2</td> </tr> <tr> <td>26-39</td> <td>± 1.5</td> </tr> <tr> <td>40-55</td> <td>± 2.4</td> </tr> </tbody> </table>	Measurement (mmHg)	Accuracy (mmHg)	5-25	± 1.2	26-39	± 1.5	40-55	± 2.4
Measurement (mmHg)	Accuracy (mmHg)								
5-25	± 1.2								
26-39	± 1.5								
40-55	± 2.4								
<p>ELECTRICAL</p> <p>Input Voltage (Tono-Pen AVIA POWERCEL Battery)</p> <p>2 x 3V Lithium Manganese Dioxide batteries (LiMnO₂)</p>	<p>OCU-FILM + TIP COVER Contains natural rubber latex</p>								
<p>SOFTWARE REVISION</p>	<p>The software revision can be obtained by contacting Reichert Technologies. The serial number identifies the manufacture date and will provide access to the software version</p>								

Guidance & Manufacturer's Declaration

Table 201 – Guidance and Manufacturer's Declaration Electromagnetic Emissions All Equipment and Systems		
Guidance and Manufacturer's Declaration – Electromagnetic Emissions		
The Tono-Pen AVIA is intended for use in the electromagnetic environment specified below. The customer or user of the Tono-Pen AVIA should ensure that it is used in such an environment.		
Emissions Test	Compliance	Electromagnetic Environment - Guidance -
RF Emissions CISPR 11	Group 1 Class B	The Tono-Pen AVIA uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
Harmonics IEC 61000-3-2	N/A	The Tono-Pen AVIA is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies building for domestic power.
Flicker IEC 61000-3-3	N/A	

Table 202 – Guidance and Manufacturer's Declaration Electromagnetic Immunity All Equipment and Systems			
Guidance and Manufacturer's Declaration – Electromagnetic Immunity			
The Tono-Pen AVIA is suitable for use in electromagnetic environment specified below. The customer or user of the Tono-Pen AVIA should ensure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
ESD IEC 61000-4-2	±6kV Contact ±8kV Air	±6kV Contact ±8kV Air	Floors should be wood, concrete or ceramic tile. If floors are synthetic, the R/H should be at least 30%.
EFT IEC 61000-4-4	±2kV Mains ±1kV I/Os	NA	Mains power quality should be that of a typical residential, commercial or hospital environment.
Surge IEC 61000-4-5	±1kV Differential ±2kV Common	NA	Mains power quality should be that of a typical residential, commercial or hospital environment.
Voltage Dips/Dropout IEC 61000-4-11	>95% Dip for 0.5 Cycle 60% Dip for 5 Cycles 30% Dip for 25 Cycles >95% Dip for 5 Seconds	NA	Mains power quality should be that of a typical residential, commercial or hospital environment. If the user of the Tono-Pen AVIA requires continued operation during power mains interruptions, it is recommended that the Tono-Pen AVIA be powered from an uninterruptible power supply or battery.
Power Frequency 50/60Hz Magnetic Field IEC 61000-4-8	3A/m	3A/m	Power frequency magnetic fields should be that of a typical residential, commercial or hospital environment.

Guidance & Manufacturer's Declaration (cont.)

Table 204 – Guidance and Manufacturer's Declaration Electromagnetic Immunity Equipment and Systems that are NOT Life-supporting			
Guidance and Manufacturer's Declaration – Electromagnetic Immunity			
The Tono-Pen AVIA is intended for use in the electromagnetic environment specified below. The customer or user of the Tono-Pen AVIA should ensure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	NA	Portable and mobile RF communications equipment should be no closer to any part of the Tono-Pen AVIA, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Radiated RF IEC 61000-4-3	80 MHz to 2.5 GHz @ 3V/m	(E1) = 3 V/m	Recommended Separation Distance: $d=(3.5/V1)(\text{Sqrt } P)$ $d=(3.5/E1)(\text{Sqrt } P)$ 80 to 800 MHz $d=(7/E1)(\text{Sqrt } P)$ 800 MHz to 2.5 GHz Where P is the max output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed transmitters, as determined by an electromagnetic site survey, should be less than the compliance levels in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol. 
Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.			
Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
* Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. The measured field strength in the location in which the ME Equipment or ME System should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the ME Equipment or ME System.			
* Over the frequency range 150 kHz to 80 MHz, field strengths should be less than $[V1] \text{ V/m}$.			

Guidance & Manufacturer’s Declaration (cont.)

Table 206 – Recommended Separation Distances between Portable and Mobile RF Communications Equipment and the Tono-Pen AVIA for ME Equipment and ME Systems that are NOT Life-supporting. Guidance and Manufacturer’s Declaration - Electromagnetic Immunity			
Recommended Separation Distances for between Portable and Mobile RF Communications Equipment and the Tono-Pen AVIA.			
The Tono-Pen AVIA is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the Tono-Pen AVIA can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF Communications Equipment and the Tono-Pen AVIA as recommended below, according to the maximum output power of the communications equipment.			
Max Output Power of Transmitter (W)	Separation (m) 150kHz to 80 MHz $d=(3.5/V1)(\text{Sqrt } P)$	Separation (m) 80 to 800 MHz $d=(3.5/E1)(\text{Sqrt } P)$	Separation (m) 800MHz to 2.5GHz $d=(7/E1)(\text{Sqrt } P)$
0.01	0.1166	0.1166	0.2333
0.1	0.3689	0.3689	0.7378
1	1.1666	1.1666	2.3333
10	3.6893	3.6893	7.3786
100	11.6666	11.6666	23.3333

Device Regulatory Classification

Insulation Protection Internally Powered (6 V battery)
 Ingress Protection IPX0
 Applied Part Type BF
 Operation Mode Continuous

Warranty

This product is warranted by Reichert, Inc. against defective material and workmanship under normal use for a period of one year from the date of invoice to the original purchaser. (An authorized dealer shall not be considered an original purchaser.) Under this warranty, Reichert's sole obligation is to repair or replace the defective part or product at Reichert's discretion.

This warranty applies to new products and does not apply to a product that has been tampered with, altered in any way, misused, damaged by accident or negligence, or which has had the serial number removed, altered or effaced. Nor shall this warranty be extended to a product installed or operated in a manner not in accordance with the applicable Reichert instruction manual, nor to a product which has been sold, serviced, installed or repaired other than by a Reichert factory, Technical Service Center, or authorized Reichert Dealer.

Lamps, bulbs, charts, cards and other expendable items are not covered by this warranty.

All claims under this warranty must be in writing and directed to the Reichert factory, Technical Service Center, or authorized instrument dealer making the original sale and must be accompanied by a copy of the purchaser's invoice.

This warranty is in lieu of all other warranties implied or expressed. All implied warranties of merchantability or fitness for a particular use are hereby disclaimed. No representative or other person is authorized to make any other obligations for Reichert. Reichert shall not be liable for any special, incidental, or consequent damages for any negligence, breach of warranty, strict liability or any other damages resulting from or relating to design, manufacture, sale, use or handling of the product.

PATENT WARRANTY

If notified promptly in writing of any action brought against the purchaser based on a claim that the instrument infringes a U.S. Patent, Reichert will defend such action at its expense and will pay costs and damages awarded in any such action, provided that Reichert shall have sole control of the defense of any such action with information and assistance (at Reichert's expense) for such defense, and of all negotiation for the settlement and compromise thereof.

PRODUCT CHANGES

Reichert reserves the right to make changes in design or to make additions to or improvements in its products without obligation to add such to products previously manufactured.

GERMANY

EICKEMEYER KG
Eltastraße 8
78532 Tuttlingen
T +49 7461 96 580 0
F +49 7461 96 580 90
info@eickemeyer.de
www.eickemeyer.de

POLAND

EICKEMEYER Sp. z o.o.
Al. Jana Pawła II 27
00-867 Warszawa
T +48 22 185 55 76
F +48 22 185 59 40
info@eickemeyer.pl
www.eickemeyer.pl

SWITZERLAND

EICKEMEYER AG
Sandgrube 29
9050 Appenzell
T +41 71 788 23 13
F +41 71 788 23 14
info@eickemeyer.ch
www.eickemeyer.ch

DENMARK

EICKEMEYER ApS
Solbakken 26, Hammelev
6500 Vojens
T +45 7020 5019
info@eickemeyer.dk
www.eickemeyer.dk

ITALY

EICKEMEYER S.R.L.
Via G. Verdi 8
65015 Montesilvano (PE)
T +39 085 935 4078
F +39 085 935 9471
info@eickemeyer.it
www.eickemeyer.it

UNITED KINGDOM

EICKEMEYER Ltd.
3 Windmill Business Village
Brooklands Close
Sunbury-on-Thames
Surrey, TW16 7DY
T +44 20 8891 2007
info@eickemeyer.co.uk
www.eickemeyer.co.uk

NETHERLANDS

EICKEMEYER B.V.
Bellweg 44
4104 BJ Culemborg
T +31 345 58 9400
info@eickemeyer.nl
www.eickemeyer.nl

CANADA

EICKEMEYER Inc.
617 Douro Street, Suite #205
Stratford, Ont. Canada
N5A 0B5
T +1 519 273 5558
F +1 519 271 7114
info@eickemeyervet.ca
www.eickemeyer.com